## Ab Maxx Pro II

Gebrauchsanleitung (Seite 2)  
Instruction manual (Page 16)  
Mode d'emploi (Page 30)  
Handleiding (Pagina 46)

DE

GB

FR

NL

Z 00246\_DE-GB-FR-NL\_V0

## Inhalt

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## Contents

Meaning of the symbols in these operating instructions	16
Safety Instructions	17
Warranty provisions	19
Product contents and product overview	20
Prior to first use	21
Unpacking the device	21
Cleaning the device	21
Inserting batteries	21
Overview of functions	22
Programme overview	23
Use	23
Preparation for using the large stimulation pad	23
Preparation for using the small stimulation pad	24
Stimulation	24
Cleaning	27
Troubleshooting	28
Technical data	29
Disposal	29
Re-ordering	29
Customer Service	29

Dear customer,

We are delighted that you have decided on the “Ab Maxx Pro II”. This special belt uses electrical muscle stimulation (EMS). Electrical pulses cause muscles to contract, which can stimulate different muscle groups at the same time. Apart from your abdominal muscles, you can also stimulate the muscles of your back, arms and shoulders.

You wish you plenty of enjoyment with your new “Ab Maxx Pro II”.

Please carefully read through the operating instructions prior to using the device for the first time and store these instructions in a safe place. The instructions are to accompany the device if it is passed on to others. The manufacturer and importer assume no liability in the event the data in these instructions has not been observed!

## Meaning of the symbols in these operating instructions

All safety notices are marked with this symbol. Please carefully read through and obey the safety notices in order to avoid injury to persons and property.

Tips and recommendations are marked with this symbol.

## Safety Instructions

### Proper Use

In order to ensure that your “Ab Maxx Pro” functions optimally, please use preferably a water-based gel as a contact medium between the silicone electrodes and your skin. You can acquire it at:

**www.service-shopping.com**

or

**+49 (0) 180-5003530** (chargeable telephone call)

- ☐ The product is intended for stimulating the muscles of your abdomen, back and shoulders.
- ☐ The product is conceived for private use, not for commercial use.
- ☐ Only use the product as described in these instructions. Any other use is deemed improper.

### Risk of Injury

- ☐ Keep pets and children away from the product and the packaging. Risk of suffocation and other dangers!
- ☐ This device is not designed to be used by persons (including children) with impaired physical, sensory or mental abilities or those lacking experience and/or knowledge unless they are supervised by a person responsible for their safety or received instructions on how to use the device from them. Children must be under supervision to ensure that they do not play with the device.
- ☐ Batteries may be life threatening if swallowed. Therefore keep batteries and the device out of the reach of young children. Should a battery be swallowed seek medical help immediately.

### Notes for Health Safety

- ☐ If you have concerns relating to your health or if you belong to any health risk group consult your physician before using the device.
- ☐ Persons with a pacemaker, heart conditions, circulatory disorders, varicose veins, diabetes, epilepsy, cancer, nerve damage or paralysis may only use the device after having consulted their physician.
- ☐ Pregnant persons should not use the device and/or only after having consulted their physician.
- ☐ The device may not be used on or by infants and sleeping or unconscious persons. You must remain awake when using the unit.
- ☐ Persons with no feeling in the skin may not use the device.
- ☐ Do not use the unit if you have a fever.
- ☐ Do not use the device if you have serious or untreated injuries (open wounds), a skin rash, pain in your calves, chronic swelling or inflamed skin.

- ☐ In the event of increasing pain or tension, turn off the device and stop using the belt. Consult a physician.
- ☐ We expressly point out that some persons may experience skin irritation or hypersensitivity brought about by the electrical muscle stimulation or the enclosed contact medium (gel). In this case, stop the use immediately and consult your physician, if necessary.
- ☐ Only use the product as indicated in the operating instructions and in the body regions stated in this manual. Any other use could result in severe damage to your health.
- ☐ Keep the product away from your head, neck, face and eyes.
- ☐ Do not use the device directly over your heart or chest.
- ☐ Do not use the product in the direct vicinity of devices that work with short or micro waves. Maintain a minimum distance of 1 meter.
- ☐ Use water-based gels as a contact medium between the silicone electrodes and your skin.
- ☐ If you experience an unpleasant prickling sensation, immediately stop using the device.
- ☐ Stimulate each muscle group no more than 30 minutes a day. We recommend limiting the use to no more than 3 times per day with at least a 4-hour break between each use. We expressly point out that improper use of the device and excessive muscle stimulation may be detrimental to your health.

#### **Danger Associated with Batteries**

- ☐ Before inserting the batteries, check whether the contacts in the device and on the battery are clean and if necessary, clean them.
- ☐ Only use the battery type specified in the technical data.
- ☐ Always replace all batteries. Do not use any different battery types, brands or batteries with different capacitance. When changing the batteries, check the polarity (+/-).
- ☐ Take the batteries out of the unit if they are empty or if you do not use the unit for a prolonged period of time. Thus you will prevent damages that could occur if the batteries leak.
- ☐ If a battery has leaked, ensure that the battery acid does not come into contact with your skin, eyes, and mucous membranes. In the event of contact with battery acid, immediately rinse the affected areas off with plenty of clean water and immediately consult a physician.
- ☐ Immediately remove a leaked battery from the device. Clean the contacts before you insert a new battery.
- ☐ The batteries may not be charged or reactivated with other means, not disassembled, thrown in a fire, submerged in a liquid or short-circuited.

#### **Device-specific safety notes**

- ☐ Only use accessories included in the product contents or exclusively those recommended by the manufacturer. Only use original spare parts provided by the manufacturer.
- ☐ Do not use the device or the belt if you experience an unpleasant prickling sensation when using it or if the belt and/or the silicone electrodes exhibit signs of heavy wear.

- ☐ If you would like to use a contact medium other than the enclosed gel, please ensure that the contact medium is preferably water-based.

#### **Environmental conditions**

- ☐ Do not subject the device to extreme temperatures, severe changes in temperatures, direct sunlight or moisture.
- ☐ Do not allow the device to fall and do not subject it to strong percussions.
- ☐ Maintain adequate distance to heat sources such as stovetops or ovens.
- ☐ Ensure that there is no risk of the device falling in water.

#### **Damages to the device**

- ☐ Check the device for damages prior to each usage session. The device may not be used if it exhibits visible damages. Do not use the device if it malfunctions or has fallen down.
- ☐ In the event the unit is defective, do not attempt to repair it yourself. In the event of damage, contact the customer service (see section "Customer service").


#### **Cleaning and storage**

- ☐ Do not use any corrosive or abrasive cleaning agents to clean the appliance. They could damage the unit.
- ☐ Only clean the housing of the device with a damp cloth. Then wipe it off completely with a dry cloth.
- ☐ Regularly clean the silicone electrodes that come into contact with your skin with a damp cloth to prevent leftover gel or residues etc. from being deposited.
- ☐ The devices and belt are not suitable for machine washing.
- ☐ Keep the device in a dry area at temperatures between 0 and 55 °C.

#### **Warranty provisions**

Defects caused by improper handling, damage or attempts at repair are excluded from the warranty. This also applies to normal wear and tear of the individual components.

## Product contents and product overview

1. Storage bag
2. Adjusting belt for the small stimulation pad (2x)
3. Small stimulation pad
4. Adjusting belt for the large stimulation pad
5. Control unit
6. Large stimulation pad

### Not illustrated:

- 1.5 volt batteries, type AAA (2x)
- Contact gel
- These instructions

## Prior to first use

Before you use the device, you must

- unpack the device and check it for transport damages,
- clean the device,
- insert the batteries in the device.

### Unpacking the device

#### ATTENTION!

- ☐ Danger of suffocation by the wrapper and bags. Keep children and pets away from the packaging.

Unpack the device and the accessories. Using the list of product contents (see “Product contents” section), check the contents for completeness and possible transport damage. Should you identify any transport damages, do not use the device (!), instead please promptly consult customer service (see section “Customer Service”).

### Cleaning the device

#### ATTENTION!

- ☐ Do not use any corrosive or abrasive cleaning agents to clean the appliance. They could damage the unit.
- ☐ Do not use any chemicals to clean the silicone electrodes.


Clean the device with a damp cloth and then immediately wipe it dry with a dry cloth.

### Inserting batteries


#### ATTENTION!

- ☐ Batteries may be life threatening if swallowed. Therefore keep batteries and the device out of the reach of young children. Should a battery be swallowed seek medical help immediately.
- You need 2 batteries for the control unit (1.5 V batteries, Type AAA). They are included in the product contents.
- When the batteries are low, the battery light will turn off. Change the batteries immediately to ensure that the device functions properly.

1. Open the battery compartment cover on the control unit by pulling it towards the arrow.
2. Insert two batteries. In the process, ensure that they are arranged with the correct polarity (+/-)!
3. Close the battery compartment.

## Overview of functions

1. With the **ON/HI** button, you turn on the device and can set the intensity of the muscle stimulation. You can choose from 10 intensity levels. The device starts at level 1. By pushing the button again, you reach the next intensity level. You have to press the button 10 times for the maximum intensity level.
2. With the **OFF/LOW** button, you can turn off the device and/or reduce the intensity. Push the button once to switch down one intensity level.
3. Select the stimulation programme with the **MODE** button. Based on the position and colour of the programme display lights, you will see which programme you have selected (see section "Programme overview").
4. The device features two integrated safety functions:
  - The device automatically deactivates if there is no contact between your skin and the silicone electrode.
  - The device automatically deactivates after 10 minutes of use.

## Programme overview

By pressing the **MODE** button, you can access the following stimulation programmes:

Programme	Name	Left programme display light	Right programme display light	Description
1	Soothing massage	RED	RED	Device alternates between 3-second muscle contractions, 2-second pause.
2	Warm-up/ warm-down	OFF	RED	Fast muscle contractions, tapping massage.
3	Muscle stimulation (normal)	GREEN	OFF	Stimulation with fast and slow components.
4	Muscle stimulation (intense)	OFF	GREEN	Alternating routine with several fast pulses.
5	Muscle stimulation (very intense)	RED	OFF	A high number of fast muscle contractions.
6	Constant stimulation	GREEN	GREEN	Constant stimulation with long muscle contractions.

## Use

### ATTENTION!

- ☐ For best results, you will need to apply a gel as a contact medium between your skin and the silicone electrodes. Make sure, however, that the contact medium is preferably water-based.
- ☐ Beginners should stimulate each muscle group for no more than 10 minutes per day!
- ☐ For beginners, we recommend limiting the use to no more than 1-2 times per day with at least a 4-hour break between each use!
- ☐ You should never exceed the maximum period of use of 30 minutes per muscle group and day. We recommend limiting use to no more than 3 workouts per day with at least a 4-hour break between workouts.

The “Ab Maxx Pro II” is not a normal massage device. It triggers muscle contractions by means of gentle electric stimulation. You can see how your own muscles move. The device does **NOT** produce any vibrations or sound waves.

### Preparation for using the large stimulation pad

1. You can use the large stimulation pad when working out your abdomen, back, shoulders, and buttocks.

2. Slide the removable large stimulation pad onto the belt.

3. Apply the control unit by using the snaps on the large stimulation pad.

4. Before you can start using the belt, you will have to apply the enclosed gel to the black silicone electrodes. The gel serves as a contact medium between the skin and electrode. Use preferably water-based gels as a contact medium.
5. Your “Ab Maxx Pro II” is now ready for use.

#### Preparation for using the small stimulation pad

1. You can use the small stimulation pad for working out muscle groups on your arm or leg.

2. Attach one of the two adjusting belts on the small stimulation pad using the Velcro® strip.
3. Secure the microprocessor on the small stimulation pad using the snap fasteners.

4. Before you can start using the belt, you will have to apply the enclosed gel to the silicone electrodes on the back of the stimulation pad. The gel serves as a contact medium between the skin and electrode. Use preferably water-based gels as a contact medium.
5. Your “Ab Maxx Pro II” is now ready for use.

#### Stimulation

1. After having realized all preparation steps, you can put on the belt.
2. Make sure that the belt is not too tight. If the belt for large stimulation pad is too small for you, you can use the elastic belts for the small stimulation pad to increase the belt circumference.
3. Arrange the silicone electrodes so that they stimulate the muscle groups that you want to stimulate. If you shift the stimulation pads during use, you will notice that muscle contractions increase or decrease. The best position is an area where you notice the strongest contractions.


4. Place the belt with the large stimulation pad on your shoulder area to relieve muscle tension by applying a soothing massage.

5. Place the belt with the large stimulation pad on your upper abdominal area in order to stimulate the muscles in this area.

6. Place the belt with the large stimulation pad on the middle of your abdominal area to stimulate the corresponding muscle groups.

7. Place the belt with the large stimulation pad on your lower abdominal area to stimulate your lower abdominal muscles.

8. Place the belt with the large stimulation pad on your waist/back area to stimulate the muscles on your side. For best results, make sure the silicone electrodes are positioned properly for this exercise! Shift the pads until you have found a position where you notice the strongest muscle contraction.

9. Place the belt with the large stimulation pad on your buttocks to firm and shape them.

10. Place the belt with the large stimulation pad on the upper or lower part of your back or in the middle of your back to stimulate your back muscles.

11. Place the belt with the small stimulation pad on the side, front, back or on the inside of your thigh to stimulate these muscle groups.

12. Place the belt with the small stimulation pad approximately 4 cm over your knee joint to stimulate the primary and knee muscles.

13. Place the belt with the small stimulation pad on your calf to stimulate the calf muscles.

14. Place the belt with the small stimulation pad on your upper arm to firm up the tissue and stimulate your biceps or triceps.

15. Place the belt with the small stimulation pad on your lower arm to stimulate your lower arm muscles.

16. After you have placed the device on the muscle group that you would like to stimulate, turn the device on with the **ON/HI** button.
17. Select a stimulation programme with the **MODE** button (see chapter “Programme overview”).
18. Select the desired intensity level with the **ON/HI** and/or **OFF/LOW** buttons.
19. The device automatically deactivates after 10 minutes.
20. Take off the device and clean it after every use as described in the section “Cleaning”.

- We recommend only using the programme 3 for your first workout week to give your muscle groups time to adjust to the new demand of the muscles.
- We recommend that beginners limit the use to a maximum of 1-2 times per day. There should be at least a 4-hour break between each use.
- When starting out, only stimulate each muscle group for no more than 10 minutes per use.
- Of course you can adjust the programme to your own level of fitness.
- We also recommend a maximum workout duration of no more than 30 minutes per muscle group and day for advanced users. Limit your use to 3 times per day with at least a 4-hour break between each use.

## Cleaning

### ATTENTION!

- ☐ Do not use any softeners, alcohol or sprays containing aerosols to clean the article.
  - ☐ Do not use any corrosive or abrasive cleaning agents to clean the appliance. They could damage the unit.
- 
1. Clean the device after every use to prevent leftover gel or other residues from being deposited.
  2. Remove the control unit from the stimulation pad before cleaning.
  3. Clean the silicone electrodes with a damp cloth and then immediately wipe them dry with a dry cloth.

## Troubleshooting

If the device is not working properly, please check whether you are able to rectify the problem yourself. Otherwise contact our customer service (see section "Customer service"). Do not attempt to repair a defective device yourself!

Problem	Possible Cause	Solution
The signal lights do not turn on.	<ul style="list-style-type: none"> <li>No batteries have been inserted.</li> <li>The batteries have not been inserted properly.</li> </ul>	<ul style="list-style-type: none"> <li>Insert batteries and press <b>ON/HI</b>.</li> <li>When inserting the batteries, ensure that the polarity is correct.</li> </ul>
The battery light turns off.	<ul style="list-style-type: none"> <li>The batteries are too weak.</li> </ul>	<ul style="list-style-type: none"> <li>Insert new batteries.</li> </ul>
The programme display lights are on, but not blinking.	<ul style="list-style-type: none"> <li>The silicone electrodes are not in contact with your skin.</li> <li>The control unit has not been properly secured to the belt.</li> <li>The intensity has been set to the lowest level.</li> <li>The control unit is not working properly.</li> </ul>	<ul style="list-style-type: none"> <li>Make sure that the stimulation pads are positioned properly and the belt is somewhat tighter.</li> <li>Check whether the control unit has been securely attached to the snap fasteners.</li> <li>Push <b>ON/HI</b> in order to increase the intensity.</li> <li>Turn the device off, take out the battery for 1 minute and attempt to restart the device.</li> </ul>
The programme display lights work, but no pulses are sent for muscle stimulation.	<ul style="list-style-type: none"> <li>The intensity level is too low.</li> <li>The belt has not been put on properly.</li> <li>The silicone electrode or stimulation pad are worn down or dirty.</li> <li>The skin is too dry.</li> </ul>	<ul style="list-style-type: none"> <li>Push the <b>ON/HI</b> button several times to increase the intensity.</li> <li>Shift the belt until you have found the position with the strongest muscle contraction.</li> <li>Replace or clean the silicone electrode or stimulation pad.</li> <li>Apply more gel on the skin.</li> </ul>

## Technical data

Model: TV90068C  
Device: Ab Maxx Pro II  
Voltage supply: 3 V DC (2x 1.5 volt battery, type AAA)  
Output: 20 ~ 180 V

## Disposal

The packaging material is recyclable. Dispose of the packaging in an environmentally-friendly manner and make it available for the recyclable material collection-service.

Dispose of the device in an environmentally-friendly manner when you decide to part with it.

The device does not belong in the household waste. Dispose of it at a recycling centre for old electric and electronic appliances.

For more information, please contact your local authority.

Notice in accordance with section 12 BattV [German Batteries Regulation]:  
Batteries and storage batteries must not be disposed of in normal household waste. As the end user you are legally responsible for returning or properly disposing of the batteries. After use you can return the batteries to the place where you bought them free of charge (e. g. at the collection point in the shop or at local authority collection points).

## Re-ordering

In order to ensure that your “Ab Maxx Pro II” functions optimally, please use preferably a water-based gel as a contact medium between the silicone electrodes and your skin. You can acquire the contact gel at our customer service or in via internet on [www.service-shopping.com](http://www.service-shopping.com)

## Customer Service

Please contact our customer service centre with any questions relating to our product:

DS Produkte GmbH, 19258 GALLIN, GERMANY

Phone: +49 (0) 180-5003530 (chargeable telephone call)

Email: [service@dspro.de](mailto:service@dspro.de)

As part of our continual expansion, we reserve the right to alter products, packaging or enclosed documentation at any time.

Importer:  
DS Produkte GmbH, Am Heisterbusch 1, 19258 GALLIN, GERMANY  
All rights reserved.

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DS Produkte GmbH, 19258 GALLIN, DUITSLAND
DS Produkte GmbH, Am Heisterbusch I, 19258 GALLIN, DUITSLAND